ASPEN® SIZED ROM KNEE

KNEE ORTHOSIS KNEIORTHESE

ÓRTESIS DE RODILLA

ORTESI DEL GINOCCHIO

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KNIE ORTHESE

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KNĚ ORTOSE

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











©2021 Aspen Medical Products, LLC. All Rights Reserved.
LB0044B revised 05/2021

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1. SIZING
A. Begin 6 inches from the middle of the patella.
B. Measure the circumference of the leg to identify the appropriate size.

SIZING GUIDE		
Size	Circumference	
X-SMALL	11 – 15 in	28 – 38 cm
SMALL	15 – 18 in	38 – 46 cm
MEDIUM	18 – 21 in	46 – 53 cm
LARGE	21 – 24 in	53 – 61 cm
X-LARGE	24 – 27 in	61 – 69 cm
XX-LARGE	27 – 35 in	69 – 89 cm

2. APPLICATION
A. To apply, open the brace and place behind the knee.
B. Wrap the upper flap around the thigh and lower flap around the calf. Ensure that the patella is centered and brace is comfortable. Be sure that the overlap is sufficient so that the hook is not exposed to the patient's skin.
C. Starting with the bottom strap, pull the hook-and-loop straps through the D-Rings and tighten to desired fit or support.
D. Secure hook-and-loop.
NOTE: Posterior straps can be adjusted for added comfort.

3. HINGE MODIFICATION
A. To adjust, remove the hinge at the top of the brace. Removal is often easier while the brace is not on the patient.
B. Use a hard surface to contour the hinge arms to maximize comfort and support.
C. To adjust the range of motion, lift the cover off the hinge.
D. Set the pins into the desired degree of flexion and extension. Replace hinge into sleeve.

INTENDED USE
To provide coronal plane stabilization of the knee while allowing sagittal plane flexion/extension to prescribed range of motion.

INDICATIONS
Knee instability due to joint laxity.

CONTRAINDICATIONS
Hypersensitivity or allergies to any of the materials from which the brace is made. Swelling of the lymphatic tissue caused by circulatory disorders. Patients who are not allowed to wear a brace according to medical instructions.

PRECAUTIONS
READ INSTRUCTIONS BEFORE USE. Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations. If you suffer from a serious medical condition, we strongly suggest that you consult with a licensed health care professional before using this product. Proper fitting is required for this product to be effective. Under some circumstances, this product may be prescribed by a physician. Please see the limited warranty for further information.

CARE INSTRUCTIONS
Hand wash only at 30° C; Wash with mild detergent; Air dry; Do not use bleach or other chemicals.

MATERIALS
Nylon, Acetal, Aluminum, Neoprene, Polyurethane Fiber (Lykra), Thermoplastic Elastomer.

NOT MADE WITH NATURAL RUBBER LATEX.

DISPOSAL
Product may be safely disposed in accordance with local laws.

REPORTING
Please report any serious incident involving the use of this device to both the Competent Authority in your state and to the manufacturer (listed on this document).

LIMITED WARRANTY
Aspen Medical Products, LLC, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products, LLC in the event of breach of warranty shall be to repair or replace the defective product or part(s).

Aspen Medical Products, LLC shall have no obligation under this limited warranty in the event: (a) The product was not purchased from Aspen Medical Products, LLC or through its authorized channels of distribution; (b) The product is altered; (c) Any parts not supplied by Aspen Medical Products, LLC are inserted into the product; or (d) The product is not used in accordance with the Aspen Medical Products, LLC Instructions for Use.

THE FOREGOING IS THE SOLE AND EXCLUSIVE REMEDY FOR FAILURE IN SERVICE OF, OR DEFECT IN, THE PRODUCT. ASPEN MEDICAL PRODUCTS, LLC SHALL NOT BE LIABLE UNDER THIS OR ANY IMPLIED WARRANTY FOR ANY DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE, AND ALL OBLIGATIONS OR LIABILITIES ON THE PART OF ASPEN MEDICAL PRODUCTS FOR DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE PRODUCTS, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED BY ASPEN MEDICAL PRODUCTS, LLC.

This warranty gives you specific legal rights and you may have other rights which vary from location to location.

A listing of authorized distributors in your area is available upon request.

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GRÖSSENTABELLE		
Größe	Umfang	
X-SMALL	11 – 15 in	28 – 38 cm
SMALL	15 – 18 in	38 – 46 cm
MEDIUM	18 – 21 in	46 – 53 cm
LARGE	21 – 24 in	53 – 61 cm
X-LARGE	24 – 27 in	61 – 69 cm
XX-LARGE	27 – 35 in	69 – 89 cm

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Dans les cas suivants, Aspen Medical Products, LLC n'aura aucune obligation au titre de cette garantie limitée:
(a) Le produit n'a pas été acheté auprès d'Aspen Medical Products, LLC ou via ses canaux de distribution autorisés;
(c) Toute pièce n'est pas fournie par Aspen Medical Products, LLC à été insérée dans le produit; ou
(d) Le produit n'est pas utilisé conformément au mode d'emploi d'Aspen Medical Products, LLC.